



Introducing the biggest surgical instrument



in orthopedic surgery.



a new level

A TABLE THAT PERFORMS LIKE A MINIMA

PROfx® ELEVATES THE SURGICAL TABLE TO A NEW LEVEL FOR

"THE PERFECT FRACTURE TABLE"

The PROfx elevates the surgical table to a new level for the orthopedic surgeon by providing the perfect table complement for any fracture procedure.

Designed to play an integral role in the surgical process, the PROfx table supports the knowledge and skill of the surgeon to minimize the invasiveness of procedures through advanced patient positioning, reliable traction capabilities and outstanding radiolucency.

ALLOWS MORE TO BE DONE THROUGH LESS INVASIVE APPROACHES.

The PROfx allows the surgeon to precisely control the position, articulation, manipulation and traction to the patient's lower extremities to an extent not possible with conventional tables. This precise maneuverability maximizes the opportunity of the chosen surgical approach to lessen the need for greater invasiveness at the incision or under the skin.

The PROfx also provides extensive imaging capability. The 35 inch cantilevered top section, combined with radiolucent leg spars provides uninterrupted imaging from chin to toes. This makes the table highly adaptable to the needs of multiple trauma patients.

With its improved patient positioning, reliable traction and outstanding radiolucency, the PROfx is ideal for orthopedic trauma cases and minimally invasive orthopedic procedures.




A REDUCTION AID IN ACETABULAR AND PELVIC FRACTURES.

Spatially guided leg spars allow unprecedented rotation, angulation and traction capabilities to enhance the surgeon's acetabular reduction options for the chosen surgical approach. Use of the PROfx table thereby limits the need for the more invasive combined or extended lateral approaches.



The table also facilitates open reduction of pelvic ring injuries and enhances the possibilities for closed reduction and percutaneous fixation.



ALLY INVASIVE SURGICAL INSTRUMENT.

OR THE ORTHOPEDIC SURGEON.

*"This type of table
has the greatest
potential for
improving results."*

- Joel M. Matta, M.D.
John C. Wilson, Jr.
Chair of Orthopedic Surgery,
Good Samaritan Hospital,
Los Angeles

res Pelvic Fractures Minimally Invasive Hip Replacement
es Hip Fractures Thoracic and Lumbar Spine Surgery



SUPPORTS MIS APPROACH TO TOTAL HIP REPLACEMENT.

The PRO_{fx} table is an essential instrument for the minimally invasive anterior approach to total hip replacement.

With its unique capability to position the leg, the PRO_{fx} enables the surgeon to replace the hip through the anterior approach without detachment of muscle from the pelvis or femur. For example, the table allows hyperextension, adduction and external rotation of the hip for femoral component placement – a positioning option not possible with conventional tables. The lack of disturbance to the lateral and posterior soft tissues provides immediate stability of the hip after surgery.



*"It is time that the
surgical table becomes
more integrated with
OR performance.
This table manages
that better than most."*

- Eric E. Johnson, M.D.
Professor of Orthopedic Surgery,
UCLA School of Medicine

IDEAL PROfx APPLICATIONS



Acetabular Fractures



Pelvic Fractures

Minimally Invasive Total Hip Replacement



Femur Fractures

Hip Fractures

Tibial Fractures

Thoracic and Lumbar Spine Surgery



PROfx RADIOLUCENCY

The combination of unrestricted C-Arm access

and a chin to toes area of unobstructed

radiolucency give the surgeon unparalleled

viewing capacity.

OR TEAM ENGINEERED



The Profx table facilitates patient positioning not

previously possible. It also streamlines OR Team

performance with several proprietary features that

make preoperative and intraoperative tasks much

easier to administer.



Features and Benefits at a Glance

Exclusive Table Maneuverability

Aids surgeon in articulation of torso and legs.
Spars include Spherical Spatial Positioning System SSPTS™.

Superior Imaging

Carbon fiber table top with extended viewing area.
Radiolucent leg spars can be positioned interoperatively for unobstructed imaging of patient's legs, hips, pelvis, lumbar and thoracic spine.

High Stability

Excellent rigidity
450 lb patient capacity
Low table profile for easy patient transfer
Automatic power floor locks enhance setup and stability.
The unique "Constant Trac" compensation mechanism built into each leg spar insures traction is maintained.

PROfx Specifications

Table top length 62 in. (158 cm)
Table length w/ spars 125 in. (316 cm)
Table top width 21 in. (53 cm)
..... 10 in. (25 cm) near perineal post
Table top height range 30 in. - 40 in. (76 cm - 102 cm)
Free and clear imaging area ... Chin to toes
Trendelenburg /
Reverse Trendelenburg 15°
Lateral tilt 20°
Leg spar articulation 35° up & down
..... 20° adduction
..... 45° abduction
Spar mount translation 8.5 in. (21 cm) each side
Patient weight capacity 450 lbs. (205 kg)

Standard Features

Carbon fiber table top
Automatic power floor locks
Standard OR accessory rail
Built-in side rails
Stainless steel base
Spherical Spatial Positioning System SSPTS™
All traction accessories for supine, prone, lateral positioning.
Safety strap
ETL approved
110V electric operating system or 220V (specify international)

Ordering Information

6850 — PROfx™ Orthopedic Surgical Table
6850I — PROfx™ Orthopedic Surgical Table (International)



Subsidiary of The Mizuho Group
30031 Ahern Avenue
Union City, CA 94587
(510) 429-1500
(800) 777-4674
Fax (888) 329-4674
Fax (510) 429-8500
www.osiosi.com

OSI is ISO certified by TUV Rheinland of North America and this table is ETL listed and CE approved.

Note: OSI is constantly improving its products, all specifications are subject to change without notice.

OSI is a Delaware Corporation. Manufactured in the USA.

©OSI NW 0428 Rev. A 3/03

